

Pressure Transducer and Block for Triaxial Test Cells

Product Code

UTGM-0110 Pressure Transducer, 2000 kPa
UTS-2405 Block with One Connection Line for Triaxial Test

The Pressure Transducer is used for the measurement of cell or back or pore pressure of water in triaxial test systems and also should be used with a U-Touch PRO Control Unit) or a UTCU-0320 (interface unit with 4 channel for data acquisition)

The Block for triaxial test cells is used for connection of the pressure transducers and de-airing in the water hoses.

