

Product Code

UTC-1070 Melting Pot, 3 L Capacity

Models for 220-240V 50-60 Hz, 1	UTC-1070
ph. Models for 110-120V 60 Hz, 1	UTC-1070-N

Standards

CEN ISO/TS 17892-2; EN 12390-3, 12390-1, 12504-1;
ASTM C31, C192, C617; AASTHO T23, T126

The Melting Pot is designed for melting capping compound, sulphur, wax and similar materials.

The melted paraffin wax is used to seal soil samples and other materials.

The apparatus consists of a 3 liter capacity aluminum container in a well-lagged steel jacket, cover and a thermostatic control heating system to keep the temperature constant in the range of ambient to 200 OC.



Dimensions	300x340x240
Weight (approx.)	mm 10 kg
Power	600 W

