

UTGE-0084 Motorized Hydraulic Specimen Extruder, Horizontal Type, 60 kN Capacity

Models for 220-240V 50-60 Hz, 1	UTGE-0084
ph. Models for 110-120V 60 Hz, 1	UTGE-0084-N

Standards

EN 13286-2, 13286-47, 12697-30; AASHTO T134, T180, T193, T245; ASTM D698, D1557, D1883, D1559; BS 1377-4, 1924-2, 598-107

The UTGE-0084 Motorized Hydraulic Extruder, Horizontal Type is designed for extruding specimens from 50 mm (2") to 125 mm (5") (outer dia.) shelly sample tubes and moulds such as proctor, CBR and marshall moulds. Also useful for sampling and removing uniaxial and triaxial test specimens to/from the cutters inner dia. from 38 mm to 100 mm.

The Extruder has 60 kN capacity and 900 mm ram travel. The hydraulic piston can be stopped at any position during the extraction.

To remove the sample from the shelly tubes, the adapter sets of different sizes on the table in next page of UTGE-0082 should be ordered separately.

The Hand Operated Hydraulic Extruder is supplied complete with

- Adapter sets for Ø100 mm (4") and Ø150 mm (6") inner dia. moulds
- Flanges and head sets for Ø38mm, Ø50mm, Ø70 mm and 100mm (only flange) inner dia. cutters.

Capacity	60kN	Weight	195 kg
Ram Travel	900 mm	(approx.) Power	750 W

Dimensions	420x2000x1200 mm
------------	------------------

For UTS-0084

UTGE-0110	Adaptor Set for 2" (50,8mm) Outer Dia Shelby Tubes.	Includes a ring with 48,5mm ID. and an extruder head with Ø43,5mm.
UTGE-0112	Adaptor Set for 2 ½" (63,5 mm) Outer Dia Shelby Tubes.	Includes a ring with 61mm ID. and an extruder head with Ø56mm.
UTGE-0114	Adaptor Set for 3" (75 mm) Outer Dia Shelby Tubes.	Includes a ring with 72mm ID. and an extruder head with Ø67mm.
UTGE-0116	Adaptor Set for 3 ¼" (83 mm) Outer Dia Shelby Tubes.	Includes a ring with 79,5mm ID. and an extruder head with Ø74,5mm.
UTGE-0118	Adaptor Set for 3 ½" (88,9 mm) Outer Dia Shelby Tubes.	Includes a ring with 85mm ID. and an extruder head with Ø80mm.
UTGE-0120	Adaptor Set for 100 mm Outer Dia Shelby Tubes.	Includes a ring with 96mm ID. and an extruder head with Ø91mm.
UTGE-0122	Adaptor Set for Ø4" (101,6 mm) Outer Dia Shelby Tubes.	Includes a ring with 97,5mm ID. and an extruder head with Ø92,5mm.
UTGE-0124	Adaptor Set for 4 ½" (114,3 mm) Outer Dia Shelby Tubes.	Includes a ring with 109mm ID. and an extruder head with Ø104mm.
UTGE-0125	Adaptor Set for 5" (125 mm) Outer Dia Shelby Tubes.	Includes a ring with 119mm ID. and an extruder head with Ø114mm.