

Product Code

UTG-0328 Motorized Sieve Shaker, 3", 5" and 8" Sieves

UTG-0329 Noise Reduction Cabinet for UTG-0323 Sieve Shaker

Models for 220-240V 50 Hz, 1	UTG-0328-T
ph. Models for 110-120V 60 Hz,	UTG-0328-N
1 ph. Models for 220-240V 60	UTG-0328-K

Standards

ASTM C136

The UTG-0328 Motorized Sieve Shaker imparts a uniform shaking action to a range of sieve dimensions. Mounting is required.



Technical Specifications

Sieve Capacity	12 pieces of 200 mm (8") sieves + pan and cover
Timer	30 minute, analogue
Power	1/4 HP

Dimensions	300x400x1150 mm
Weight (approx.)	30 kg

