

Product Code

UTG-0180 Ultrasonic Cleansing Apparatus 28

L,
220-240 V 50-60 Hz, 1 ph.

UTG-0185 Ultrasonic Cleansing Apparatus 50

L, 220-240 V 50-60 Hz, 1 ph.

The UTG Series Ultrasonic Cleansing Apparatus is used for efficient and safe cleaning of sieves suitable for fine mesh sieves which could be damaged using ordinary cleaning methods. The tank and tank components are manufactured from stainless chrome.

Time and temperature can be easily adjusted by the digital control unit.



Technical Specifications

UTGG-0180	
Capacity	28 L
Ultrasonic Power	1200 W Peak / 600 W Effective
Heater Power	1200 W
Dimensions (External)	330x530x500 mm
Dimensions (Internal)	200x300x500 mm
Weight (approx.)	10 kg

UTGG-0185	
Capacity	50 L
Ultrasonic Power	2000 W Peak / 1000 W
Heater Power	Effective 1200 W
Dimensions (External)	600x650x980 mm
Dimensions (Internal)	300x400x500 mm
Weight (approx.)	20 kg