

Product Code

UTW-0900 Calibration Weight Set M1 Class, (1g-1kg) S. Steel
UTW-0901 Calibration Weight, 500 g, M1 Class, Cast Iron
UTW-0902 Calibration Weight, 1 kg, M1 Class, Cast Iron
UTW-0903 Calibration Weight, 2 kg, M1 Class, Cast Iron
UTW-0904 Calibration Weight, 5 kg, M1 Class, Cast Iron
UTW-0905 Calibration Weight, 10 kg, M1 Class, Cast Iron
UTW-0906 Calibration Weight, 20 kg, M1 Class, Cast Iron
UTW-0907 Calibration Weight, 25 kg, M1 Class, Cast Iron
UTW-0920 Calibration Weight, Set F1 Class, (1g-1kg) S. Steel

The UTEST range of calibration weights are used for the periodical verification of your laboratory balances.

The UTW-0900 M1 Class and UTW-0920 F1 Class Calibration Weight Sets are manufactured from stainless. Except weight set, other M1 class weights are manufactured from cast iron coated with black epoxy based paint.

Traceable calibration certificates should be ordered separately if required.

The Calibration Weight Sets are supplied complete with a wooden box.



Calibration Weight Sets (UTW-0900 M1 Class and UTW-0920 F1 Class)

Weigh	1 g	2 g	5 g	10 g	20 g	50 g	100 g	200 g	500 g	1 kg
t Unit	1	2	1	1	2	1	1	2	1	1