

Product Code

**UTW-0800 Mechanical Balance 310 g x 0.01 g**

**UTW-0810 Mechanical Balance 2610 g x 0.1 g**

The UTW Series Mechanical Balances, with their metal base construction, beam design and stainless steel pans, are durable, precise and easy-to-use.

They are especially designed for on-site applications when electric power is not available. All models are equipped with zero adjust knobs at the end of the beam for fast and easy zeroing.

The UTW-0800 model eliminates the need for level adjustment with its special three-point base design to perform high-precision measurements.

The UTW-0810 is known for its reliable performance at a reasonable price.



UTW-0800

	UTW-0800	UTW-0810
Capacity	310 g	2610 g
Readability	0.01 g	0.1 g
Platform Size	13x89 mm	147x147 mm
Dimensions (mm)	200x550x220	200x550x220
Weight (approx)	5 kg	5 kg



UTW-0810