

Product Code

- UTW-0633 Digital Balance 600 g x 0.01 g**
- UTW-0635 Digital Balance 6 kg x 0.1 g**
- UTW-0637 Digital Balance 30 kg x 5 g**
- UTW-0638 Digital Balance 30 kg x 0.5 g**
- UTW-0642 Digital Balance 3000 g x 0.01 g**
- UTW-0643 Digital Balance 4500 g x 0.01 g**
- UTW-0644 Digital Balance 6200 g x 0.01 g**
- UTW-0645 Digital Balance 15 kg x 0.2 g**
- UTW-0648 Digital Balance 30 kg x 0.1 g**
- UTW-0654 Digital Balance 60 kg x 1 g**

The UTW Series Digital Balances provide a wide range of maximum capacity and readability characteristics which make them economical and easy to use, they are ideal for central and site laboratories who require a range of balances for various applications.

The UTW Series Digital Balances are fitted with strain gauge load cells and are designed with large backlit LCD displays that give precise measurements within the 0°C to 40°C temperature range. All models can be connected to printers or PC's through their RS 232 outputs.



UTW-0635

	Capacity	Readability	Pan Size	Additional Power Supply
UTW-0633	600 g	0.01 g	Ø122 mm	DC 6V (300mA) Batt. (AA)
UTW-0635	6 kg	0.1 g	322x212mm	220-240V
UTW-0637	30 kg	5 g	226x187mm	Pb Batt. and Adapter
UTW-0638	30 kg	0.5 g	322x212mm	220-240V
UTW-0642	3000 g	0.01 g	Ø150mm	DC 12V and Battary (AA)
UTW-0644	6200 g	0.01 g	195x195mm	DC 12V (1,2A)
UTW-0645	5 kg	0.2 g	306x222mm	Pb Batt. and Adapter
UTW-0648	30 kg	0.1 g	322x212mm	220-240V
UTW-0654	60 kg	10 g	500x600mm	Pb Batt. and Adapter

	Under-Bench Weighing Facility	Dimensions	Weight (approx.)
UTW-0633	No	210x310x140mm	2 kg
UTW-0635	Yes	400x400x150mm	4,50 kg
UTW-0637	No	400x400x150mm	4,50 kg
UTW-0638	Yes	400x400x150mm	4,50 kg
UTW-0642	Yes	370x470x340mm	5 kg
UTW-0643	Yes	370x470x340mm	5 kg
UTW-0644	No	400x400x150mm	4,50 kg
UTW-0645	Yes	400x400x150mm	4,50 kg
UTW-0648	No	400x400x150 mm	4,50 kg