

Product Code

UTW-0622 Electronic Analytical Balance, LCD Screen
210 g x 0.0001 g

UTW-0625 Electronic Analytical Balance 200 g x 0.001 g

UTW-0628 Electronic Analytical Balance 300 g x 0.001 g

UTW-0630 Electronic Analytical Balance 500 g x 0.001 g

UTW-0631 Electronic Analytical Balance 820 g x 0.001 g

UTW-0632 Electronic Analytical Balance 1000 g x 0.001 g

The UTW Series Electronic Analytical Balances set the standard in terms of reliable measurements, user-friendly and long service life. The balances can be used for routine tasks and/or complicated weighing procedures.

The UTW-0622, UTW-0625 and UTW-0628 models have an automatic internal calibration feature, the other models in the range require external calibration.

The UTW Series Analytical Balances provide fast and easy installation and cleaning with their ergonomic designs and brilliant displays. They can also be connected to printers or PC's through RS 232 outputs to give excellent monitoring and reporting performance.



UTW-0622

	Capacity	Readability	Pan Size	Additional Power Supply
UTW-0622	220 g	0.1 mg	Ø90 mm	DC 12V (1,2A) Adapter
UTW-0625	210 g	1 mg	Ø115 mm	DC 12V (1,2A) Adapter
UTW-0628	300 g	1 mg	Ø115 mm	DC 12V (1,2A) Adapter
UTW-0630	500 g	1 mg	Ø115 mm	DC 12V (1,2A) Adapter
UTW-0631	820 g	1 mg	Ø115 mm	DC 12V (1,2A) Adapter
UTW-0632	1000 g	1 mg	Ø115 mm	DC 12V (1,2A) Adapter

	Under-Bench Weighing Facility	Dimensions	Weight (approx.)
UTW-0622	No	400x500x500 mm	6 kg
UTW-0625	No	250x250x300 mm	2 kg
UTW-0628	Yes	250x250x300 mm	2 kg
UTW-0630	Yes	250x250x300 mm	2 kg
UTW-0631	Yes	380x470x340 mm	5 kg
UTW-0632	Yes	380x470x340 mm	5 kg