

Product Code

UTD-1415 Warm-Air Drier, 220-240 V 50-60 Hz

UTD-1418 Hot Air Gun, 220-240 V 50-60 Hz

The UTD-1415 Warm-Air Drier and UTD-1418 Hot-Air Gun are used for drying small amounts of aggregate particles and soil samples.

Both models have the heat and air flow control option. The UTD-1418 has 300-500°C air temperature and 240-450 L/min air flow capacity.

	UTD-1415	UTD-1418
Dimensions	240x280x100 mm	200x300x100mm
Weight (approx.)	1 kg	0,5705 kg W
Power	2200 W	1600 W



UTD-1418



UTD-1415