

## Product Code

**UTB-1500 Apparatus for Distillation of Cutback Asphalt**  
**UTGT-2025 High Distillation Thermometer,**  
**Range-2 +400°C, ASTM 8C**

Models for 220-240V 50-60 Hz, 1 ph.	UTB-1500
Models for 110-120V 60 Hz, 1 ph.	UTB-1500-N

## Standards

**ASTM D402; AASHTO**  
**T78**

The UTB-1500 is used for the examination of cut-back asphaltic materials by the distillation test. The apparatus consists of a distillation flask, condenser, adapter, shield, electric heater with thermoregulator, cylinder receiver, thermometer (ASTM 8C), shield and flask support.

Dimensions	300x300x600 mm
Weight (approx.)	6 kg

