

Product Code

UTB-1350 Cleveland Flash Tester, 220-240 V 50-60 Hz
UTGT-2040 Thermometer IP28C, -6 +400°C

Standards

EN 2592; ASTM D92; AASHTO T48

The UTB-1350 Cleveland Flash Tester is used for determining the flash and fire point of petroleum products.

It consists of a brass cup mounted on an electric heater with a temperature controller and a thermometer (-6 C to +400 C).

Conforming to the CE European Directive.



Dimensions	380x380x380 mm
Weight (approx.)	4 kg
Power	600 W