

Product Code

UTS-0160 Soil Lathe / Trimmer and Extruder UTS-0164
Wire Saw with one Blade
UTS-0165 Wire Saw Blade for UTS-0164 (pack of 10)
UTS-0166 Trimming Knife
UTGG-2205 Porcelain Mortar with Pestle 130 mm dia
UTGG-2215 Rubber Headed Pestle

The UTS-0160 Soil Lathe / Trimmer and Extruder is used to extrude and trim soil samples from 35 mm to 100 mm diameter to reduce samples.

Wire Saw, Trimming Knife, Porcelain Mortar with Pestle, the Rubber Headed Pestle can be ordered separately.



UTS-0160 with UTS-0164 and UTS-0166

Technical Specifications

Specimen Lathe	35x70 mm to 100x200
Specimen Trimming and Extrusion	mm 35x70 mm to 50x100
Vertical Daylight	mm 260 mm

Dimensions	230x300x430
Weight (approx.)	mm 14 kg