

Product Code

UTGE-4300 pH Indicator Papers, pH Range 1 to 14
 UTGE-4320 Quantab Chloride Titrator Type 1175, 40 strips/pack
 UTGE-4322 Quantab Chloride Titrator Type 1176, 40 strips/pack

Standards

BS 812:117, 1377:3

The UTGE-4300 pH Indicator Papers are used for quick determination of pH in the 1 to 14 pH range.

UTGE-4320 and UTGE-4322 Quantab Chloride Titrators are used for quick determination of water soluble chloride salts present in soils and aggregates. It is based on the Volhard Method.

UTGE-4320 covers 0.005% to 0.1% NaCl and UTGE-4322 covers 0.05% to 1% NaCl.



UTGE-4300

Product Code	Dimensions	Weight (approx.)
UTGE-4300	50x10x80 mm	0.1
UTGE-4320	70x70x120	kg
UTGE-4322	mm	0.1
	70x70x120	kg
	mm	0.1
		kg



UTGE-4320 & UTGE-4322