

Product Code

UTS-0280 Ion Exchange Apparatus
UTS-0282 Ion Exchange Resin 500 g

Standards

BS 1377:3

The UTS-0280 Ion Exchange Apparatus when used together with UTS-0282 Ion Exchange Resin, is used to determine the sulphate content of aqueous soil extracts and ground water. The apparatus consists of an ion exchange column of 10 mm diameter and 400 mm long, swan-neck outlet and a 1500 ml round bottom flask to give a constant head. The apparatus is supplied assembled on a stand.

UTS-0282 Ion Exchange Resin, 500 g should be ordered separately.

	UTS-0280	UTS-0282
Dimensions	200x100x600 mm	50x50x50 mm
Weight (approx.)	5 kg	0,50 kg

