

Product Code

UTCM-0003 Le Chatelier Flask

Standards

EN 196-6, 450-1, 15617-1; ASTM C110, C128, C188; C989; AASHTO T133

The UTCM-0003 Le Chatelier Flask is used to determine the density of hydraulic cement, ground granulated blast-furnace slag and fly ash for concrete, filler aggregates and lime. The glass flask has a 250ml capacity. The neck is graduated from 0 to 1 ml and from 18 to 24 mL in 0.1-mL graduations.



Dimensions	100x100x300
Weight (approx.)	mm 0.1 kg