

### Product Code

UTCM-0010 Le Chatelier Mould  
UTCM-0014 Le Chatelier Soundness Kit


### Standards

EN 196-3, 450-1, 459-2; EN ISO 9597



UTCM-0010

The soundness of cement, fly ash for concrete and lime is determined by using the Le Chatelier moulds and Le Chatelier Water Bath (UTCM-0016) according to the relevant standard.



UTCM-0014

The Le Chatelier Soundness Kit is supplied complete with

- Le Chatelier moulds 3 pcs.
- 50x50 mm glass plates 6 pcs.
- 300 gr Weights 1 pcs.
- 100 gr Weight, 3 pcs.
- Tamping Rod 17 mm dia. x 70 gr
- Steel Ruler
- Plastic Carrying Case

Dimensions	340x290x80
Weight (approx.)	mm 2 kg