

Product Code

UTCM-0020 High Pressure Cement Autoclave,
230 V 50-60 Hz, 1ph

Standards

ASTM C151, C490; UNE 7207

The UTCM-0020 High Pressure Cement Autoclave is designed to perform expansion tests on cement specimens. 10 specimens can simultaneously be tested in the high pressure steam vessel of 154 mm diameter and 430 mm height.

The Autoclave consists of a pressure gauge, pressure regulator, temperature regulator, control switches safety valve and specimen rack. Certified conforming to INAIL procedure.

UTCM-0033A Two Gang Shrinkage mould and Length Comparators (UTCM-0200 or UTCM-0210) should be ordered separately.

The High Pressure Cement Autoclave is supplied complete with

- Specimen Rack, 10 samples capacity



Dimensions	450x475x1080
Weight (approx.)	mm 55 kg
Power	2600 W