

Product Code

UTCM-0875E Manuel Mortar Mixer, EN
 UTCM-0875A Manuel Mortar Mixer, ASTM
 UTCM-0876 Spare Bowl, 5 L
 UTCM-0878E Spare Beater, EN
 UTCM-0878A Spare Beater, ASTM

Models for 220-240V 50-60 Hz, 1 ph.	UTCM-0875E	UTCM-0875E
	UTCM-0875E-N	UTCM-0875E-N

Standards

EN 196-1, 196-3, 413-2, 459-2, 480-1, 1015-2, 12617-4;
 ASTM C187, C305, AASHTO T129, T131, T162

The UTCM-0875 Manual Mortar Mixer has a 5 liter (approx.) capacity, they have been designed to mix mortars and cement pastes primarily to the requirements of standards.

The mixing paddles revolve at a rate of 140 rpm. with a planetary motion of 62 rpm. in low speed. In high speed, the paddles revolve at the rate of 285 rpm. with a planetary motion of 125 rpm.

The user can choose the speeds easily by using the switch fitted to the machine. There is a sand filling apparatus on the mixers to pour sand easily. The bowl and beater are easily fitted and removed from mixer.

The Manual Mortar Mixers are supplied complete with; •
 • Bowl, 5 L (approx.)
 • Beater



Dimensions	300x555x610
------------	-------------

Weight (approx.)	mm 54 kg
------------------	----------

Power	200 W
-------	-------



UTCM-0878



UTCM-0876