

Product Code

UTAS-0195 Asphalt Mixer, 7.5 L
UTAS-0196 Spare Mixing Bowl, 7.5 L, for UTAS-0195
UTAS-0197 Spare Mixing Whisk, for UTAS-0195

Standards

EN 12697-35

The mixer has a capacity of 7.5 liters and is designed for mixing of asphalt samples to be used for mechanical tests such as compaction, indirect tensile, and Marshall. The mixer does not contain a heater.

The Laboratory Mixer is supplied complete with

- Bowl, 7,5 lt Capacity, 5,5 kg, Stainless Steel
- Mixing Whisk

Dimensions	370x570x700
Weight (approx.)	mm 70 kg
Power	550 W



UTAS-0196



UTAS-0197

