

Product Code

- UTGT-1430 Viscosity Thermometer +19 to +27°C 0.1°C Div., ASTM 17C
- UTGT-1432 Viscosity Thermometer +34 to +42°C 0.1°C Div., ASTM 18C
- UTGT-1434 Viscosity Thermometer +49 to +57°C 0.1°C Div., ASTM 19C
- UTGT-1436 Viscosity Thermometer +57 to +65°C 0.1°C Div., ASTM 20C
- UTGT-1438 Viscosity Thermometer +79 to +87°C 0.1°C Div., ASTM 21C
- UTGT-1440 Viscosity Thermometer +95 to +103°C 0.1°C Div., ASTM 22C
- UTGT-1500 Mechanical Hygrometer/Thermometer/Barometer
- UTGT-1520 Mechanical Hygrometer
- UTGT-1550 Mechanical Dial Thermometer 0 to +260°C
- UTGT-1580 Digital Stop Watch

UTEST Viscosity Thermometers are used for UTAS-0300 Saybolt Two-Tube Digital Viscometer in order to determine the Saybolt Viscosity of petroleum products at specified temperatures where accurate temperature measurements are required.

UTGT-1500 Mechanical Hygrometer/Thermometer/Barometer is a reliable and useful instrument to monitor the humidity, temperature and the atmospheric pressure simultaneously.

UTGT-1520 Mechanical Hygrometer is reliable instrument for measuring the humidity of the ambient air and UTGT- 1550 Mechanical Dial Thermometer provides accurate temperature measurements especially at high temperatures up to 260°C.

UTGT-1580 Digital Stop Watch is designed for accurate time measurements with 0.01 second sensitivity.



Viscosity Thermometers



UTGT-1520



UTGT-1550



UTGT-1500



UTGT-1580

Product Code	Description	ASTM Ref.	IP Ref.	Dimensions	Weight
UTGT-1430	Viscosity Thermometer +19 to +27°C 0.1°C Div.	17 C	-	30x30x300 mm	0.1 kg
UTGT-1432	Viscosity Thermometer +34 to +42°C 0.1°C Div.	18 C	23 C	30x30x300 mm	0.1 kg
UTGT-1434	Viscosity Thermometer +49 to +57°C 0.1°C Div.	19 C	-	30x30x300 mm	0.1 kg
UTGT-1436	Viscosity Thermometer +57 to +65°C 0.1°C Div.	20 C	-	30x30x300 mm	0.1 kg
UTGT-1438	Viscosity Thermometer +79 to +87°C 0.1°C Div.	21 C	-	30x30x300 mm	0.1 kg
UTGT-1440	Viscosity Thermometer +95 to +103°C 0.1°C Div.	22 C	24 C	30x30x300 mm	0.1 kg
UTGT-1500	Mechanical Hygrometer/Thermometer/Barometer	-	-	150x150x150 mm	0.2 kg
UTGT-1520	Mechanical Hygrometer	-	-	150x150x150 mm	0.2 kg
UTGT-1550	Mechanical Dial Thermometer 0 to +260°C	-	-	50x50x200 mm	0.2 kg
UTGT-1580	Digital Stop Watch	-	-	75x65x25 mm	0.1 kg