

Product Code

- UTGE-3700 Laboratory Air Compressor 8 bar - 25 L
- UTGE-3705 Laboratory Air Compressor 8 bar - 50 L
- UTGE-3710 Laboratory Air Compressor 10 bar - 100 L

Models for 220-240V 50-60 Hz, 1 ph

UTGE-3700	UTGE-3705	UTGE-3710
-----------	-----------	-----------

UTGE series Laboratory Air Compressors are durable machines for supplying compressed air required by several analyses in material testing.

The compressors are supplied complete with;

- Air Gun
- Hose

	Voltage	Dimensions	Weight (approx)
UTGE-3700	220-240 V 50-60 Hz	600x300x600 mm	28.5 kg
UTGE-3705	220-240 V 50-60 Hz	600x300x600 mm	40 kg
UTGE-3710	220-240 V 50-60 Hz	1460x470x860 mm	140 kg

*110 V 60 Hz should be mentioned on order if required

	Motor Power	Tank Capacity	Air Suction Capacity	Max. Working Pressure
UTGE-3700	2 HP	25 L	192 L/min	8 bar
UTGE-3705	2 HP	50 L	192 L/min	8 bar
UTGE-3710	2 HP	100 L	250 L/min	10 bar



Air Gun and Hose



UTGE-3700



UTGE-3710