

Product Code

UTG-0326 Rotatap Motorized Sieve Shaker, with Time Adjustment, for 200 mm (8") Dia. Frame Sieves

UTG-0327 Noise Reduction Cabinet for UTG-0326 Sieve Shaker

Models for 220-240V 50 Hz, 1 ph.	UTG-0326-T
Models for 110-120V 60 Hz, 1 ph.	UTG-0326-N
Models for 220-240V 60 Hz, 1 ph.	UTG-0326-K

Standards

ASTM C136

UTG-0326 Rotatap Motorized Sieve Shaker is designed to be used for 8" diameter sieves. This shaker provides 278 oscillations with a horizontal circular motion and 150 taps per minute with vertical tapping motion to produce an effective sieving action.



Technical Specifications

Sieve Capacity	7 pieces of 200 mm (8") sieves with full height + pan 14 pieces of 200 mm (8") sieves with full height + pan
Timer	99 minute, digital
Power	1/4 HP
Dimensions	510x700x640 mm
Weight (approx.)	75 kg